

Fluke Biomedical Ansur

Test and Inspection Procedure

Copyright © 2000 - 2012 Fluke Biomedical

Test record

TEST PASSED

Test performed

Date: 19/9/2022
Record: IEC 62353 Alternative
Leakage - CL1.mtr
Template: IEC 62353 Alternative
Leakage - CL1.mtt

Ansur components used

Ansur Version 2.9.7
Plug-In: ESA620 Version 1.1.11

Test setup

Selections

Service events performed

Standards performed

IEC 62353 (CL1)

Device under test

Serial number	24084	Type	Centrífuga
Appliance code		Model	DCA-300RVT
Group	Laboratorio	Location	Laboratorio
Status	MP	Address 1	Catamarca
Manufacturer	Presvac	Address 2	2do Piso

MTI Data

Test instrument	Serial number	Firmware version
ESA 620	2982052	v2.10

Signatures

Test result

Test element		Test type				Fail
IEC 62353 Alternative Leakage - Class I		Auto Sequence				
Mains Voltage		Mains Voltage				
Live to Neutral	Result:	Value	Unit	High limit	Low limit	Standard
	Live to Neutral	221,3	V			IEC 62353 (CL1)
Neutral to Earth	Result:	Value	Unit	High limit	Low limit	Standard
	Neutral to Earth	222,8	V			IEC 62353 (CL1)
Live to Earth	Result:	Value	Unit	High limit	Low limit	Standard
	Live to Earth	0,8	V			IEC 62353 (CL1)
Protective Earth Resistance		Protective Earth Resistance				
Configuration:						
Test Current: Low						
Alternative Equipment Leakage		Alternative Equipment Leakage				
Closed Earth	Result:	Value	Unit	High limit	Low limit	Standard
	Closed Earth	676,0	uAAC	1000		IEC 62353 (CL1)
Open Earth	Result:	Value	Unit	High limit	Low limit	Standard
	Open Earth	676,0	uAAC	1000		IEC 62353 (CL1)
Insulation Resistance		Insulation Resistance				
Configuration:						
Test Voltage: 500V						
Mains to Protective Earth	Result:	Value	Unit	High limit	Low limit	Standard
	Mains to Protective Earth	99999,0	MOhm			IEC 62353 (CL1)